

ORCA-Flash4.0 V2 Digital CMOS Camera C11440-22CU/C11440-22CU01 Instruction manual

Thank you for your purchase

 CAUTION	<ul style="list-style-type: none">• Follow the safety precautions in Chapter 1 in order to avoid personal injury and damage to property when using this camera. Be sure to read this Instruction manual beforehand in order to use the digital camera correctly. The manual describes the correct method of handling the camera and provides cautions in order to avoid accidents.• After reading, keep the manual where it can be referred to at any time.
--	--



Ver.1.7
Sept 2016

HAMAMATSU PHOTONICS K.K.

1. SAFETY PRECAUTIONS







1-1 INDICATION OF THE SYMBOLS

The following symbols can be found on this camera:

	Direct current
	Alternating current

1-2 CLASSIFICATION OF WARNING

We have classified the warnings symbols that appear in this instruction manual and on the camera as follows for your convenience. Make sure that you fully understand them and obey the instructions they contain.

 WARNING	Improper handling of the camera without observing these warnings could lead to serious injury to the user and even death.
 CAUTION	Improper handling of the camera without observing these cautions could lead to personal injury to the user or damage to property.
 Note	This symbol indicates a note to help you get the best performance from the camera. Read the contents of the note carefully to ensure correct and safe use. Failure to observe one of these notes might impair the performance of the camera.
	This symbol indicates a cautionary item that should be obeyed when handling the camera. Read the contents carefully to ensure correct and safe use.
	This symbol indicates an action that is forbidden. Read the contents carefully and be sure to obey them.
	This symbol indicates a compulsory action or instruction. Read the contents carefully and be sure to obey them.

 **WARNING**



Power supply

Use the camera with the voltage indicated on the rating sticker. Using a different voltage can damage the camera and lead to fire or electric shock.



Cables

Be careful not to place heavy objects on cables or bend it excessively. Doing so can damage the cable and lead to fire or electric shock.



Power supply cord

Use the accessory power supply cord when this camera is used.



AC adapter

Use the accessory AC adapter when this camera is used.



Do not touch the plug with wet hand. Doing so can lead to electric shock.



Do not attempt to dismantle or modify the camera

Doing so can also lead to damage and even injury, as some internal components become very hot. Only touch parts as indicated in this manual.



Do not insert a foreign substance into the camera

Do not allow foreign objects such as combustible substances, metal objects or water to get inside the camera. They can damage the camera and lead to fire or electric shock.



If an abnormality occurs

Such as the image suddenly disappearing or a strange noise, smell or see smoke coming from the camera, stop the power supply immediately and contact Hamamatsu subsidiary or local distributor. Never attempt to repair the camera yourself.

 **CAUTION**



AC adapter

When unplugging the AC adapter, do not pull on the cord, but remove the plug from the outlet to avoid causing electric shock or fire.



When unplugging the AC adapter, do not pull on the cord, but remove the plug from the camera to avoid breakdown of the AC adapter or the camera.



When the camera is not used for a long period of time, unplug the power supply cord from the outlet to avoid damaging the cord and causing electric shock or fire.



Connecting and disconnecting cables

Always turn off the power supply of the peripheral device before connecting and disconnecting cables.



Fixed the camera

When fitting the camera to a tripod or other fixture, use the optional base plate. Be careful that the fitting screw does not enter more than 8 mm from the surface of the base plate. Screwing it in excessively can impair normal operation.



Lenses (C11440-22CU)

Be careful not to screw the lens more than 7 mm onto the C-mount of the camera. Doing so can scratch the protective glass. (Some wide-angle lenses in particular can have a thread of 7 mm or more.)



Shipping precautions

When transporting the camera by truck, ship, airplane, etc., wrap it securely in packaging material or something similar.



Strong impact

Do not subject the camera to strong shocks by dropping it, for example. Doing so can damage the camera.



At the water cooling

Be careful water does not splash on the camera. Cut off the power supply of the circulating water cooler and the camera when you remove and install the cooling water hoses.



Operating environment

This camera is designed and tested for use in an industrial environment. If this camera is used in residential areas, EMI (electro-magnetic interference) may occur. This camera must not be used in residential areas.

⚠ CAUTION



Disposal

When disposing of the camera, take appropriate measures in compliance with applicable regulations regarding waste disposal and correctly dispose of it yourself, or entrust disposal to a licensed industrial waste disposal company. In any case, be sure to comply with the regulations in your country, state, region or province to ensure the camera is disposed of legally and correctly.




Cooling water

It is recommended to use soft water (except pure water) for cooling water. Follow instruction manual which is attached to your circulating water cooler for an appropriate temperature range of cooling water. If you plan on using water other than soft water as recommended for example antifreeze etc, please refer to description of cooling water which is written in 12. [Maintenance] or contact Hamamatsu subsidiary or local distributor.



Condensation

At the Water-cooling, if ambient temperature and ambient humidity become high, condensation will take place easily. Use the camera under the environment where condensation will not take place referring to the following graph.



2. CHECK THE CONTENTS OF PACKAGE

When you open the package, check that the following items are included before use. If the contents are incorrect, insufficient or damaged in any way, contact your local dealer without attempting to operate the camera.

Camera: C11440-22CU or C11440-22CU01	1
AC adapter	1
Power supply cord for AC adapter	1
Lens mount cap (attached to the camera)	1
C11440-22CU/C11440-22CU01 Before Use (Booklet)	1
C11440-22CU/C11440-22CU01 Instruction manual (CD-ROM)	1
QC sheet	1

[Option]

Cooling water hose (2 hoses)	A10788-04
SMA-BNC cable	A12106-05
SMA-SMA cable	A12107-05
Camera Link interface cable	A11255-05
USB 3.0 interface cable	A12467-03
USB 3.0 interface board	M9982-25
Adjuster pole for C11440-22C	A11185-01
Base plate for C11440-22C	A11186-01
Camera Link interface board	M9982-20

Note

- The cable listed in option is highly recommended for use with the camera. The camera system may not confirm to CE marking regulation if other type of cable is used with.

Note

- The cooling is set as air cooling at factory setting. It is possible to change the cooling setting from Air cooling to Water cooling by the special software included.



- Handle the circulating water cooler and the cooling water according to an instruction manual of the circulating water cooler.

Note

- If you use the adjuster pole and the base plate, see each installation manual.

3. INSTALLATION

Avoid using or storing this camera in the following places



- Where can be out of the specification of the operating temperature.
- Where can be out of the specification of the storage temperature.
- Where the temperature varies extremely.
- In direct sunlight or near a heater.
- Where can be out of specification of the operating humidity and can be dripped with liquid.
- Where can be out of specification of the storage humidity and can be dripped with liquid.
- Close to a strong source of magnetism or radio waves.
- Where there is vibration.
- Where it might come into contact with corrosive gases. (such as chlorine or fluorine)
- Where there is a lot of dust.


How to place the camera (when the camera is placed on a table)



Place the camera the water connectors to be lateral side.



Do not place the camera the rear panel of the camera, which connectors are located, to be at the bottom. (Do not block ventilation openings.)



Do not block ventilation openings



To prevent overheating in the camera's interior, do not wrap the camera in cloth or other material, or in any way allow the camera's ventilation ports to become blocked. If the camera is being operated in an enclosed environment, ensure clearance of at least 2 cm from both the intake and exhaust vents when setting up.

Weight of the camera



Be careful not to drop off the camera or not drop underfoot when making it move because it is approx. 2 kg.

Contents

1. SAFETY PRECAUTIONS	1
1-1 INDICATION OF THE SYMBOLS.....	1
1-2 CLASSIFICATION OF WARNING.....	1
2. CHECK THE CONTENTS OF PACKAGE	5
3. INSTALLATION	6
4. OVERVIEW	9
5. FEATURES	9
6. NAME AND FUNCTION OF THE PARTS	12
7. CONNECTION	15
8. WATER COOLING	18
8-1 CAUTIONS.....	18
8-2 CONNECTION OF WATER COOLING HOSES.....	20
8-3 DISCONNECTION OF WATER COOLING HOSES.....	21
9. OPERATION	22
9-1 PRECAUTIONS.....	22
9-2 PREPARATION FOR IMAGING.....	23
9-2-1 WHEN USING AIR-COOLING.....	23
9-2-2 WHEN USING WATER-COOLING.....	23
9-3 END OF IMAGING.....	23
10. DESCRIPTION OF VARIOUS FUNCTIONS	24
10-1 THEORY OF CMOS IMAGE SENSOR.....	24
10-2 NORMAL AREA MODE.....	26
10-2-1 READOUT METHOD (SCAN MODE).....	26
10-2-2 FRAME RATE CALCULATION.....	27
10-2-3 CONFIGURING EXPOSURE TIME.....	31
10-2-4 FREE RUNNING MODE.....	32
10-2-5 EXTERNAL TRIGGER MODE.....	34
10-2-6 TRIGGER OUTPUT.....	39
10-2-7 SLOW SCAN.....	41
10-2-8 GLOBAL RESET.....	42
10-3 LIGHT SHEET READOUT MODE.....	44
10-3-1 READ-OUT DIRECTION.....	44
10-3-2 FRAME RATE CALCULATION.....	46
10-3-3 READ-OUT TIME OF THE HORIZONTAL LINE.....	46
10-3-4 TIMING DIAGRAM.....	47
10-3-5 PROGRAMMABLE TIMING OUTPUT.....	49
10-4 REAL-TIME CORRECTION FUNCTIONS.....	50

11. PRECAUTION WHEN USING CMOS IMAGE SENSOR.....	51
12. MAINTENANCE.....	52
12-1 CARE.....	52
12-2 INFORMATION ON COOLING WATER FOR THE CIRCULATING WATER COOLER.....	53
12-2-1 WHEN USING COOLING WATER OTHER THAN RECOMMENDED.....	53
13. TROUBLESHOOTING CHECKLIST	54
13-1 IMAGE IS NOT TRANSFERRED.....	54
13-2 ALTHOUGH IMAGES ARE TRANSFFERED	54
14. SPECIFICATIONS.....	55
14-1 CAMERA SPECIFICATIONS	55
14-2 CONDENSATION.....	58
14-3 SPECTRAL RESPONSE CHARACTERISTICS(TYP.).....	59
14-4 INTERFACE SPECIFICATIONS.....	60
14-4-1 CAMERA LINK INTERFACE	60
14-4-2 OUTPUT TIMING SPECIFICATIONS	63
15. DIMENSIONAL OUTLINES.....	67
15-1 C11440-22CU.....	67
15-2 C11440-22CU01.....	68
16. WARRANTY	69
17. CONTACT INFORMATION	70

4. OVERVIEW

C11440-22CU/C11440-22CU01 is equipped with the new scientific image sensor, an advanced CMOS device that realizes the multiple benefits of high resolution, high readout speed, and low noise all at once.

The camera provides 4.0 megapixels resolution at 100 frame/s (and up to 25 655 frame/s by sub-array readout) while achieving 1.0 electrons (median) 1.6 electrons (r.m.s) readout noise performance. Moreover, the camera delivers high sensitivity through its on-chip micro lens, 37 000:1 high dynamic range that make the camera suitable for almost any scientific application from bright field imaging to low-light fluorescence imaging across a wide spectral range. Various external trigger functions and timing output functions ensure proper timing control with peripheral equipment to cover a wide range of applications.

The camera is the new scientific digital camera for life science microscopy, semiconductor inspection, x-ray scintillator readout or industrial imaging.

5. FEATURES

(1) Readout noise

In the camera, the pixel amplifier is optimized: it has high gain from optimizing the semiconductor process, and the difference among pixel amplifiers are greatly minimized. In addition, there is on-chip CDS (correlated double sampling) circuit, which plays an important role in achieving low noise. Moreover, the sensor features a split readout scheme in which the top and bottom halves of the sensor are readout independently, and the data of each horizontal line is read by 2 lines of column amplifier and A/D in the top and the bottom in parallel and simultaneously. As a result, it achieves very fast readout speed while keeping very good low-noise performance.

The camera has lower readout noise (1.0 electrons (median), 1.6 electrons (r.m.s)) than the conventional cooled CCD camera. Moreover, high-speed readout (100 frame/s with 2048 pixels × 2048 pixels) with very low readout noise, which was impossible, can now be achieved.

(2) Cooling structure

In the camera, CMOS image sensor is cooled down by the peltier element to suppress the dark current. If CMOS image sensor is exposed to the atmosphere, condensation of the moisture from the air might occur. However the camera has a special hermetic chamber structure to isolate the sensor from the atmosphere, and the chamber is filled with nitrogen gas.

(3) Pixel number and pixel size

CMOS image sensor has 6.5 μm × 6.5 μm pixel sizes that is equivalent to conventional CCD image sensor (2/3 inch, 1.3 megapixels). Also, the camera can observe a wider field of view because the pixel number is about 3 times that of the conventional CCD image sensor (2/3 inch, 1.3 megapixels)

(4) Readout method

The camera has a variety of readout modes. In addition to full resolution readout mode, sub-array readout and binning readout are supported.

(5) Frame rate (readout speed)

CMOS image sensor which this camera adopts realizes both low noise (1.0 electrons (median) 1.6 electrons (r.m.s)) and high speed readout (100 frame/s with 2048 pixels x 2048 pixels) simultaneously, by a split readout scheme in which the top and the bottom halves of the sensor are readout independently, and the data of each horizontal line is read by 2 lines of column amplifier and A/D in the top and the bottom in parallel and simultaneously.

(6) Real-time correction functions

When using the camera, there is a case that shading caused by uneven illumination or optics is not negligible in the image. Also, there are a few pixels in CMOS image sensor that have slightly higher readout noise performance compared to surrounding pixels. For those cases, the camera has real-time offset level, shading and defective pixel correction features to further improve image quality. The correction is performed in real-time without sacrificing the readout speed at all.

(7) Interface

This camera has both Camera Link and USB 3.0 interface.

■ Camera Link Interface

Camera Link interface is able to transfer large volumes of data. It can transfer 4 megapixels image with 100 frame/s.

In order to realize such large volume data transfer, the camera uses Camera Link "Full Configuration Deca Mode" which is an expanded version of Camera Link "Full Configuration". It enables to transfer maximum 85 MHz x 10 Taps (8 bit) image data to computer as fast as 100 frame/s.

In order to use this interface, a Camera Link interface board which supports "Full Configuration Deca Mode" is required.

■ USB 3.0 Interface

USB 3.0 interface is able to transfer 4 megapixels image with 30 frame/s. It is versatile interface which is suitable to use when fast data transfer is not required.

This interface does not require a Camera Link interface board. It transfers image with moderate transfer speed.




• Do not connect Camera Link and USB 3.0 interface simultaneously.



• When a connection interface is changed from Camera Link to USB 3.0, and vice versa, the application software must be closed and the camera must be turned off.


(8) Camera readout modes

The camera has two scan speed in Normal area mode, and two readout direction in Light sheet readout mode.




(9) Camera operation modes

The camera has three operation modes: 1) the free running mode, in which the exposure and readout timing are controlled by the internal microprocessor, and 2) the external trigger mode, in which the exposure and readout timing are decided by an external trigger. 3) the start trigger mode is used to start operating the camera by a trigger input for a continuous imaging.





6. NAME AND FUNCTION OF THE PARTS

(1) C11440-22CU



(2) C11440-22CU01



 CAUTION	<ul style="list-style-type: none"> Place the camera the water connectors to be lateral side. Do not place the rear panel of the camera, which connectors are located, to be at the bottom (Do not block ventilation openings.)
--	---


① **Lens mount**

C11440-22CU can be attached to C-mount lens or an optics system.
 C11440-22CU01 can be attached to F-mount lens or an optics system.

Note	<ul style="list-style-type: none"> The depth of the C-mount is 7 mm. Screwing in the mount too far can scratch the glass surface.
-------------	--

② **WATER connector [WATER] (at Water-cooling)**

It connects the camera and the circulating water cooler with the cooling water hoses. The insert position of WATER IN/OUT is not specified.

 CAUTION	<ul style="list-style-type: none"> See 8. [Water cooling] for details.
--	---

③ **Camera Link interface connector [DIGITAL OUT 1]**

This is connected to the Camera Link interface connector 1 on the computer.

④ **Camera Link interface connector [DIGITAL OUT 2]**

This is connected to the Camera Link interface connector 2 on the computer.

⑤ **Trigger input connector [EXT.TRIG]**

This is used when the camera is being operated using external synchronization.
 Input is 3.3 V LVCMOS level, and input impedance is 10 kΩ.
 When an external trigger is input, the trigger is activated at the falling or rising edge of the signal. (You can choose external trigger polarity between Negative and Positive.)

⑥ **Timing out connector 1,2,3 [TIMING 1,2,3]**


This is used when peripheral device(s) require synchronization with the camera.
 Output is 3.3 V LVCMOS level, and it is output though BUS TRANSCEIVER IC SN74AVC8T245.
 Output impedance is 33 Ω.

Note	<ul style="list-style-type: none"> Determine termination according to cable length and so on.
-------------	--

⑦ **STATUS lamp [STATUS]**

The LED indicates status of camera.

Lighting color	Status of power distribution
Turn off (no color)	Power off
Green (Blinking)	Initialization
Green (lighting)	Power on
Orange (lighting)	Data transfer
Red (lighting)	Heat up

 CAUTION	<ul style="list-style-type: none"> When the camera heats up, stop operation and unplug the AC adapter immediately.
--	---

⑧ **Power switch [POWER]**

The power is turned on/off.

When the power switch is set to "ON", the camera turns on and starts initialization and the lamp blinks in green.

When the initialization is completed, the lamp color stays in green.

When the camera transfers data and the lamp color is orange.

When the power switch is set to "OFF", the camera returns to the power off state and the lamp turns off.

⑨ **DC power input connector [DC IN]**

This is the power supply terminal. Use the accessory AC adapter.

⑩ **USB 3.0 interface connector [USB 3.0]**

This is connected to the USB 3.0 interface connector on the computer.



- Do not connect Camera Link and USB 3.0 interface simultaneously.



- When a connection interface is changed from Camera Link to USB 3.0, and vice versa, the application software must be closed and the camera must be turned off.

⑪ **Air inlet**

This is the inlet for the heat ventilation.



- To prevent overheating inside the camera, do not wrap the camera in cloth or other material, or block the camera's ventilation.



- If the camera is being operated in an enclosed environment, ensure to keep clearance at least 2 cm from both intake and exhaust vents when setting up.

⑫ **Screw holes for attaching the option**

They are used for attaching an adjuster pole(A11185-01) or a base plate(A11186-01). A11185-01 and A11186-01 are the options of C11440-22C.



- Refer to each manual about the attaching method of the option.

7. CONNECTION



- When you connect cables, turn off the power supply of the camera and the peripheral devices.

Note

- If you use the adjuster pole and the base plate, see each installation manual.

CAUTION

- Place the camera the water connectors to be lateral side. Do not place the rear panel of the camera, which connectors are located, to be at the bottom (Do not block ventilation openings.)

Refer to the figure when connecting the various cables.

■ Overall connection




Figure 7-1

■ Connection to Camera Link




Figure 7-2

■ Connection to USB 3.0




Figure 7-3



- Do not connect Camera Link and USB 3.0 interface simultaneously.



- When a connection interface is changed from Camera Link to USB 3.0, and vice versa, the application software must be closed and the camera must be turned off.

① AC adapter

This is the cord to supply a power supply. Use the accessory AC adapter.

② Cooling water hose (at Water-cooling: Option)

It connects the camera and circulating water cooler. The insert position of WATER IN/OUT on the camera WATER connector is not specified.



CAUTION

- See 8. [Water cooling] for details.

③ **Camera Link interface cable 1 (Option)**

This is the cable to connect the Camera Link interface connector 1 of the camera and the Camera Link interface connector 1 on the computer.

④ **Camera Link interface cable 2 (Option)**

This is the cable to connect the Camera Link interface connector 2 of the camera and the Camera Link interface connector 2 on the computer.

Note

- Hamamatsu recommends A11255-05 optional Camera Link interface cable for this camera. The camera complies with EMC direction with using A11255-05 Camera Link interface cable. Be careful that the camera with other interface cable may not fulfill the EMC directive requirements.

⑤ **USB 3.0 interface cable (Option)**

This is the cable to connect the USB 3.0 interface connector of the camera and the USB 3.0 interface connector on the computer.

Note

- Hamamatsu recommends A12467-03 optional USB 3.0 interface cable for this camera. The camera complies with EMC direction with using A12467-03 Camera Link interface cable. Be careful that the camera with other interface cable may not fulfill the EMC directive requirements.

8. WATER COOLING



CAUTION

- Improper handling of the camera without observing these cautions could lead to personal injury to the user or damage to property.

8-1 CAUTIONS

(1) Cooling water


It is recommended to use soft water (except pure water) for cooling water. If you plan on using water other than soft water as recommended for example antifreeze etc, please refer to description of cooling water which is written in 12. [Maintenance] or contact Hamamatsu subsidiary or local distributor.

(2) Recommendation ambient temperature

Hamamatsu recommends 20 °C for Circulating water temperature. For the appropriate temperature range of the cooling water, confirm with the instruction manual of your circulating water cooler.

(3) Condensation

Use the camera under the environment where condensation will not take place referring to the following graph.



(4) Handling of the circulating water cooler

Handle the circulating water cooler and the cooling water according to an instruction manual of the circulating water cooler.

Proper performance may not be achievable if a non-recommended circulating water cooler is used.

(5) Start water cooling and water cooling in operation

- Confirm the water is flowing before starting the camera cooling and that the camera does cool.
- Keep 0.45 L/min flow rate for water circulation.
- Do not stop coolant while the camera is working.

(6) Cooling water hose

The hose has a blue hose (Internal diameter: 4 mm / External diameter: 6 mm) and a gray hose (Internal diameter: 8 mm / External diameter: 13.5 mm). (Figure 8-1)

If the hose size on circulating water cooler is the same as blue hose, remove gray hose from the joint part. The gray hose can be removed when blue hose is pulled with pushing the button of the joint on gray hose. (Figure 8-2)




Figure 8-1




Figure 8-2

(7) Connection of the cooling water hose

CAUTION

- Follow the instruction in Section 8-2 [CONNECTION OF WATER COOLING HOSES] and Section 8-3 [DISCONNECTION OF WATER COOLING HOSES] to connect / disconnect the hose.

- Stop water circulation when connecting / disconnecting the hose, and turn off the power of the camera and the circulating water cooler.
- Confirm that cooling water stops.
- Prepare water absorption sheet (such as Waste, Towel or so) and catch pan in order to avoid water drop or water splash.

(8) Deterioration of the cooling water hose

Replace the water hose with a new one whenever it cannot keep 0.45 L/min flow rate for water circulation due to the hose deterioration.

8-2 CONNECTION OF WATER COOLING HOSES




Figure 8-3

- (1) Place the camera on the stable table.
- (2) Connect water cooling hose into the WATER connector on the camera.
 - Insert the hose fully into the WATER connector on the camera. (as shown in Figure 8-3)
 - Confirm the hose stops at it.
- (3) Set the camera onto a microscope (If the camera is used on the microscope).

If it is easy to connect the hose onto the camera after the camera is set onto the microscope then it is OK to connect the hose after the camera is set on the microscope.
- (4) Connect the hose onto the circulating water cooler.

Follow the instruction on the circulating water cooler when you connect the hose onto the circulating water cooler.
- (5) Turn on the circulating water cooler and confirm the cooling water is flowing normally.

CAUTION

- Stop the circulating water cooler when the water flow is abnormal or water drop or splash is found.

8-3 DISCONNECTION OF WATER COOLING HOSES



- Remove the water cooling hoses only when it is necessary to remove.

- Turn off the camera power and all peripheral devices including circulating water cooler.
- Remove the hose on circulating water cooler side.
Follow the instruction on the circulating water cooler when you disconnect the hose from the circulating water cooler.
- Remove water or water drop inside the hose and the camera by air.
 - Blow Air from one side of hose. Prepare water absorption sheet (such as Waste, Towel or so) and catch pan on another side of hose in order to avoid water drop or water splash.
 - Blow Air until no water drop come out.
- Remove the camera from the microscope (if the camera is used on the microscope).
It is not necessary to remove the camera from the microscope if it is possible to remove the hoses from the camera as it is.
- Place the camera on the stable table.
Put the lens cap on to protect the sensor.
- Change the WATER connector direction to be downward.
Prepare water absorption sheet (such as Waste, Towel or so) and catch pan.
- Remove hoses one by one, and wipe water.
Disconnect hoses with pushing button while being careful not to splash water.



Figure 8-4

CAUTION

- Cooling water may be left inside the camera even after hoses are removed. In such case, remove water inside by blowing air from connectors. Be careful not to splash water onto the camera.

9. OPERATION

9-1 PRECAUTIONS

Be careful of the following when you operate the camera.

(1) Cooling method

Cooling of this equipment is done using a Peltier element. With a Peltier element, when current is supplied, one surface is cooled, and the other surface is heated. CMOS image sensor is positioned on the cooled side, and cooling is done by discharging the heat from the heated surface.

The camera has two cooling modes, Air-cooling mode and Water-cooling mode. Cooling mode can be changed by software which is called, "DCAM Configurator". After cooling mode was changed, the camera memorizes the last setting as the default setting for cooling. The present cooling mode set-up of this camera can be checked using "DCAM Configurator".

Note

- A software "DCAM Configurator " is in a CD-ROM of DCAM software.

Cooling method	Detail
Air-cooling mode (Forced air-cooled) (Default)	The heated side of a peltier element is cooled by a fan inside the camera. When the camera is turned on, the fan starts rotating and cooling is started.
Water-cooling mode	Circulating water cooler (Optional) is used for cooling the heated side of a peltier element. Cooling does not start just turning on the camera. Cooling water circulation must be started before start operating the camera in water-cooling. A fan inside the camera does not rotate. (Please refer to 8 [WATER COOLING] for instruction of water-cooling.)



- Do not switch to water-cooling mode when water-cooling is unnecessary.

(2) Ambient temperature

The recommended ambient temperature for camera operation is 20 °C. Both water-cooling or air-cooling are available as cooling method, CMOS image sensor cooling temperature is more stable under water cooling operation.

(3) Protection circuit

This camera's thermoelectric cooling device is protected by a protection circuit. If the internal of the camera becomes abnormally hot, the protection circuit operates to inform that by a buzzer alarm (beep tone) and lighting the red LED light while simultaneously cutting the current supply to the Peltier element. As soon as this protection is implemented, turn off the power switch. Then remove the cause of the overheating..

9-2 PREPARATION FOR IMAGING

Use the following procedure when start operating the camera.



- When you connect cables, turn off the power supply of the camera and the peripheral devices.



- Change cooling method (cooling mode) if needed. The present cooling mode set-up of this camera can be checked using DCAM Configurator.

9-2-1 WHEN USING AIR-COOLING

- (1) Connect the equipment as shown in Figure 7-1 before operating of the camera.
- (2) Turn on the camera.
- (3) Check cooling fan is operating properly and air is circulating.

9-2-2 WHEN USING WATER-COOLING

- (1) Connect the equipment as shown in Figure 7-1 before operating of the camera.
- (2) Turn on the camera.
- (3) Turn on the circulating water cooler.
- (4) Check cooling water is circulating suitably.
- (5) Turn on the cooling switch of the camera from application software.

Note

- Refer to the manual of application software for ON/OFF of the cooling switch of a camera.

9-3 END OF IMAGING

Follow the procedure below when imaging is finished.

- (1) End the imaging or transmission of image data with the application software.
- (2) Turn off the circulating water cooler. (In the case of water-cooling)
- (3) Turn off the camera and peripheral devices.

10. DESCRIPTION OF VARIOUS FUNCTIONS

10-1 THEORY OF CMOS IMAGE SENSOR

The pixel of a CMOS image sensor is composed of the photodiode and the amplifier that converts the charge into voltage. Entered light is converted to charge and converted to voltage in the pixel. The voltage of each pixel is output by switching the switch one by one. (Figure 10-1)

The scientific CMOS image sensor used in this camera has an on-chip CDS (correlated double sampling) circuit, which plays an important role in achieving low noise. In addition, CMOS image sensor realizes both low noise and high speed readout simultaneously, by a split readout scheme in which the top and the bottom halves of the sensor are readout independently, and the data of each horizontal line is read by 2 lines of column amplifier and A/D in the top and the bottom in parallel and simultaneously.




Figure 10-1 Structure of CMOS image sensor

The exposure and the readout method of CMOS image sensor which this camera adopts is rolling shutter. In the rolling shutter, the exposure and readout are done line by line. Therefore, the exposure timing is different on one screen. (Figure 10-2) But even if the object moves during the exposure, the affect of rolling shutter is very small.




Figure 10-2 Readout timing of Rolling shutter

10-2 NORMAL AREA MODE

10-2-1 READOUT METHOD (SCAN MODE)

(1) Readout direction

The camera readout the image sensor from the center line to the top and from the center line to the bottom simultaneously (center line is depicted in red line in the diagram).




Figure 10-3 Normal area mode readout direction

(2) Normal readout

Perform charge readout from camera individually for all pixels.

(3) Binning readout

With this camera, 2x2 binning readout and 4x4 binning are available by adding the signal of adjacent pixels in the digital domain, Binning readout is a method for achieving high sensitivity in exchange for losing resolution.

(4) Sub-array readout

Sub-array readout is a procedure only a region of interest is scanned. It is possible to increase the frame rate by reducing the number of vertical lines scanned. When a target area is placed in the center of the screen, sub-array readout can perform the fastest readout. In sub-array readout, binning configuration is enabled.

Size and a position of the readout area can be configured according to the table below.

Interface	Binning	Settings			
		Horizontal		Vertical	
		Size	Position	Size	Position
Camera Link	Regardless of binning	---	---	4 lines steps	
USB 3.0	1x1 (Normal readout)	512 pixels	32 pixels	8 lines	4 lines
	2x2 binning readout	256 pixels	16 pixels	4 lines	2 lines
	4x4 binning readout	128 pixels	8 pixels	2 lines	1 lines

Note

- Refer to 10-2-2 [FRAME RATE CALCULATION] about the frame rate of each readout mode.

10-2-2 FRAME RATE CALCULATION

10-2-2-1 At Standard Scan

(1) Camera Link

V_n = Number of vertical line (The center of the set area is the middle of the sensor.)
 $Exp1$ = 1.004 ms to 10 s
 $1H$ = $9.744\ 36 \times 10^{-6}$



- The Exp1 value must be input to the calculation formula below in units of seconds.

[Calculation formula, Value of Frame rate]

	Calculation formula	Horizontal / Vertical		Frame rate (frame/s)
		Horizontal	Vertical	
Free running mode	$1/(V_n/2 \times 1H)$	2048	2048	100
			1024	200
			512	400
			256	801
			128	1603
			64	3206
			8	25 655
External trigger mode (Edge trigger / Level trigger)	$1/(V_n/2 \times 1H + Exp1 + 1H \times 10)$	2048	2048	90
			1024	164
			512	278
			256	425
			128	579
			64	707
			8	876
External trigger mode (Synchronous trigger readout)	$1/(V_n/2 \times 1H + 1H \times 18)$	2048	2048	98
			1024	193
			512	374
			256	702
			128	1251
			64	2052
			8	4664
External trigger mode (Fast Synchronous trigger readout)	$1/(V_n/2 \times 1H + 1H \times 5)$	2048	2048	99
			1024	198
			512	393
			256	771
			128	1487
			64	2773
			8	11 402

(2) USB 3.0

- V_n = Number of vertical line (The center of the set area is the middle of the sensor.)
 $Exp1$ = 1.004 ms to 10 s
 $Exp2$ = 1.05 ms to 10 s
 $1H$ = $9.744\ 36 \times 10^{-6}$
 $int()$ = The decimal point is rounded down.



- The $Exp1$ and $Exp2$ values must be input to the calculation formula below in units of seconds.

[Calculation formula]

	Binning	Horizontal \ Vertical		Calculation formula
		Horizontal	Vertical	
Free running mode	1x1	1024	---	$1 / (int(V_n/2048/30/1H) \times 1H)$
		1536		
		2048		
	1x1	512	---	$1 / (int(V_n/2048/100/1H) \times 1H)$
	2x2	---	---	
4x4	---	---		
External trigger mode (Edge trigger)	1x1	1024	$112 \leq V_n \leq 2048$	$1 / (int(V_n/2048/30/1H) \times 1H)$
		1536		
		2048		
	1x1	1024	$8 \leq V_n \leq 104$	$1 / (V_n/2 \times 1H + 1H \times 10 + Exp1)$
		1536		
		2048		
1x1	512	---	$1 / (V_n/2 \times 1H + 1H \times 10 + Exp1)$	
2x2	---	---		
4x4	---	---		
External trigger mode (Level trigger)	1x1	1024	$112 \leq V_n \leq 2048$	$1 / (int(V_n/2048/30/1H) \times 1H)$
		1536		
		2048		
	1x1	1024	$8 \leq V_n \leq 104$	$1 / (V_n/2 \times 1H + 1H \times 4 + Exp2)$
		1536		
		2048		
1x1	512	---	$1 / (V_n/2 \times 1H + 1H \times 4 + Exp2)$	
2x2	---	---		
4x4	---	---		
External trigger mode (Synchronous trigger readout)	1x1	1024	$24 \leq V_n \leq 2048$	$1 / (int(V_n/2048/30/1H) \times 1H)$
		1536		
		2048		
	1x1	64	$8 \leq V_n \leq 16$	$1 / (V_n/2 \times 1H + 1H \times 18)$
		1024		
		1536		
1x1	2048	---	$1 / (V_n/2 \times 1H + 1H \times 18)$	
1x1	512	---		
2x2	---	---		
4x4	---	---		

[Value of frame rate (unit: frame/s)]

	Binning	1×1 (Normal readout)		2×2 / 4×4
	Horizontal Vertical	1024 / 1536 / 2048	512	---
Free running mode	2048	30	100	100
	1024	60	200	200
	512	120	400	400
	256	240	801	801
	128	481	1603	1603
	64	968	3206	3206
	8	7894	25 655	25 655
External trigger mode (Edge trigger)	2048	30	90	90
	1024	60	164	164
	512	120	278	278
	256	240	425	425
	128	481	579	579
	64	707	707	707
	8	876	876	876
External trigger mode (Level trigger)	2048	30	90	90
	1024	60	164	164
	512	120	279	279
	256	240	428	428
	128	481	583	583
	64	713	713	713
	8	886	886	886
External trigger mode (Synchronous trigger readout)	2048	30	98	98
	1024	60	193	193
	512	120	374	374
	256	240	702	702
	128	481	1251	1251
	64	968	2052	2052
	8	4664	4664	4664

Note

- The calculation formula and the frame rate value of Start trigger mode (External trigger mode) are same as the free running mode. About this mode, refer to 10-2-5-4 [Start trigger mode].

10-2-2-2 At Slow Scan

(1) Camera Link, USB 3.0

- V_n = Number of vertical line (The center of the set area is the middle of the sensor.)
 $Exp1$ = 3.021 ms to 10 s
 $1H$ = 32.4812×10^{-6}



- The Exp1 value must be input to the calculation formula below in units of seconds.

[Calculation formula, Value of Frame rate]

	Calculation formula	Horizontal / Vertical		Frame rate (frame/s)
		Horizontal	Vertical	
Free running mode	$1/(V_n/2 \times 1H)$	2048	2048	30
			1024	60
			512	120
			256	240
			128	481
			64	962
			8	7696
External trigger mode (Edge trigger / Level trigger)	$1/(V_n/2 \times 1H + Exp1 + 1H \times 10)$	2048	2048	27
			1024	50
			512	85
			256	133
			128	184
			64	228
			8	287
External trigger mode (Synchronous trigger readout)	$1/(V_n/2 \times 1H + 1H \times 18)$	2048	2048	29
			1024	58
			512	112
			256	210
			128	375
			64	615
			8	1399
External trigger mode (Fast Synchronous trigger readout) (Camera Link)	$1/(V_n/2 \times 1H + 1H \times 5)$	2048	2048	29
			1024	59
			512	117
			256	231
			128	446
			64	832
			8	3420

10-2-3 CONFIGURING EXPOSURE TIME

The exposure time setting can be done by absolute value. The actual exposure time setting is defined by the following formula, and the camera automatically calculates a larger and closest value from the specified exposure time setting.



- The Exp1 values must be input to the calculation formula below in units of seconds.

(1) Standard scan

Exp1 = 1 ms to 10 s (at sub-array 38.96 μ s to 10 s)

Exp2 = $\text{Exp1} \times 10^6 \div 9.744\ 36\ \mu\text{s}$ (round up at decimal point)

$$9.744\ 36\ \mu\text{s} \times \text{Exp2}$$

(2) Slow scan

Exp1 = 3 ms to 10 s (at sub-array 129.99 μ s to 10 s)

Exp2 = $\text{Exp1} \times 10^6 \div 32.4812\ \mu\text{s}$ (round up at decimal point)

$$32.4812\ \mu\text{s} \times \text{Exp2}$$

Available setting range of the exposure time is the following.

	Standard scan	Slow scan
Free running mode	1 ms to 10 s	3 ms to 10 s
Free running mode (at Sub-array)	38.96 μ s* to 10 s	129.99 μ s to 10 s
External trigger mode	1 ms to 10 s	3 ms to 10 s

Note

* 38.96 μ s (Standard scan) and 129.99 μ s (Slow scan) is the minimum exposure time when sub-array is set to 8 lines vertically symmetric (4 lines in top half and 4 lines in bottom half) with respect to the horizontally center axis. The minimum exposure time vary depend on vertical line number of sub-array setting.

10-2-4 FREE RUNNING MODE

The camera has the free running mode which the exposure and readout timing can be set by software command and controlled by an internal microprocessor. The free running mode has normal readout mode (in which the exposure time is longer than the 1 frame readout time) and electrical shutter mode (in which the exposure time is shorter than the 1 frame readout time). These readout modes are automatically switched depending on the exposure time setting.

Note

- Contact to Hamamatsu subsidiary or local distributor for the detail of the timing information.

10-2-4-1 Normal readout mode

The normal readout mode is suitable for observation, monitoring, field of view and focus adjustment, and animation because it can operate with full resolution, which is faster than the video rate*. (* at Camera Link = 100 frame/s, at USB 3.0 = 30 frame/s)

In addition, the exposure time can be extended to collect more signals and increase the signal-to-noise ratio if the object is dark. In the normal readout mode, the exposure time is the same or longer than the 1 frame readout time. In this mode, the frame rate depends on the exposure time, and it becomes frame rate = 1/exposure time. The maximum exposure time is 10 s.




Figure 10-4

10-2-4-2 Electrical shutter mode

The electrical shutter mode is used to get a proper signal level when signal overflow happens due to too much input photons in normal readout mode. In this mode, the fastest frame rate is 100 Hz (via Camera Link) or 30 Hz (via USB 3.0) at full resolution even when the exposure time is short.

[Camera Link]




Figure 10-5

[USB 3.0]




Figure 10-6

10-2-5 EXTERNAL TRIGGER MODE

The camera has various external trigger functions to synchronize the camera with the external equipment. In the external trigger mode, the external equipment becomes a master and the camera becomes a slave.

Note

- Contact to Hamamatsu subsidiary or local distributor for the detail of the timing information.

10-2-5-1 Edge trigger mode

The edge trigger mode is used so that the exposure starts according to an external signal. Exposure time is set by software command. In this mode, the exposure of the first line begins on the edge (rising/falling) timing of the input trigger signal into the camera. (1023H and 1024H in the following figure) The exposure of the second line is begun after the readout time of one line passes (1022H and 1025H in the following figure), and the exposure is begun one by one for each line.




Figure 10-7 (Ex. rising edge)

10-2-5-2 Level trigger mode

The level trigger mode is used to control both exposure start timing and exposure time length by inputting external trigger pulses. In the mode, the camera starts exposure at the start of high or low period of the input trigger pulse and stops exposure at the end of high or low period of the input trigger pulse. The example below is for the trigger level High. The exposure of the first line begins when the trigger signal becomes High, and the exposure of the second line begins after the readout time of line one passes. Each exposure begins one by one for each line. The exposure of the first line is finished when the trigger signal becomes low, and signal readout is begun. The exposure time of each line is defined by the time that the input trigger is high. The minimum trigger pulse width is 1 ms + 50 μ s.




Figure 10-8 (Ex. rising edge)

10-2-5-3 Synchronous readout trigger mode

The synchronous readout trigger mode is used for continuous imaging when it is necessary to control the exposure start timing of each frame from an external source. It is useful for confocal microscopy. For example, when the camera is used with a spinning disk confocal microscope and the camera exposure time is synchronized to the spinning disk's rotation speed, it is possible to eliminate uneven illumination (called banding noise) caused by variation of the spinning disk rotation speed. Also, it is useful for securing as long exposure time as possible while controlling the exposure start timings by external trigger signals.

It is possible to work with shorter interval time if you set fast synchronous readout trigger mode. Fast synchronous readout trigger mode works with only Camera Link interface operation.

(1) Normal operation (when the pulse count is set as 1.)

The synchronous readout trigger mode is used for continuous imaging when it is necessary to control the exposure start timing of each frame from an outside source and also when it is necessary to secure as long exposure time as possible. In the synchronous readout trigger mode, the camera ends each exposure, starts the readout and also, at the same time, starts the next exposure at the edge of the input trigger signal (rising / falling edge). That is, the interval between the same edges of the input trigger becomes the exposure time.




Figure 10-9 (Ex. rising edge)

(2) Pulse count

Also in the synchronous readout trigger mode, synchronous readout can be controlled by specifying, set by command, the number of timing pulses to determine the exposure time. The following figure shows the exposure timing when the pulse count is set as 3.




Figure 10-10 (Pulse count)

10-2-5-4 Start trigger mode

The start trigger mode is to start operating the camera by a trigger input for a continuous imaging. It is useful to secure the frame rate as fast as possible when continuous image acquisition and not to sacrifice the exposure time. For example, when it is necessary to measure the phenomenon after stimulation, it is possible to start continuous image acquisition at the stimulation timing.

The start trigger mode is to start operating the camera by a trigger input for continuous imaging, and it works at the highest frame rate because it is operated in free running mode. In the start trigger mode, the camera starts exposure and switches to free running mode by the edge of an external trigger signal (rising / falling edge).




Figure 10-11 (Ex. rising edge)

10-2-5-5 External trigger delay function

In most cases when a delay is needed between the laser pulse emission and the exposure start is needed, a delay unit is set between the laser and camera to control trigger timing. In each external trigger mode of the camera, the delay can be set to the trigger signal input to the camera by command. With this setting, a range of trigger can be arranged without a delay unit. The range for delay time is 0 μ s to 10 s (10 μ s steps).

10-2-6 TRIGGER OUTPUT

The camera provides a range of trigger output signals to synchronize with an external instrument and the camera becomes the master and the external instrument becomes the slave. There are three different trigger output functions as follows. Also, it can output continuous High output (High output fixed) or continuous Low output (Low output fixed).

These three different trigger output functions can be selected by software command, and they are output from Timing out connector.

Note

- Refer to Figure 10-4 to Figure 10-11 about details of each trigger output functions.

10-2-6-1 Global exposure timing output

It shows the global exposure timing where all lines expose at the same time. There is a case that one event is divided into two frames because the timing of the exposure in each line is different for the rolling shutter. However, by using the Global exposure timing output the global exposure becomes possible for the phenomenon that happens for this period. Global exposure timing output shows the period where all lines expose at the same time.

Note

- There is no output signal when the exposure time is less than the frame rate.

10-2-6-2 Programmable timing output

By using the programmable timing output, synchronizing external devices is simple. A system that needs simple timing signal does not require a delay unit or pulse generator. It is possible to program and output a pulse that has an optional pulse width and an optional delay time to the end of readout timing or Vsync. The setting range for delay time is 0 μ s to 10 s, and the setting range for pulse width is 1 μ s to 10 s.

The relation between the parameter which can be set with each reference signal, and an output signal becomes below.

Reference signal	Output signal
Read End	Camera outputs a pulse after certain delay, from the end of sensor readout. Also the pulse width can be set.
Vsync	Camera outputs a pulse after certain delay, from the beginning of readout. Also the pulse width can be set.




Figure 10-12 Programmable timing output

10-2-6-3 Trigger ready output

The trigger ready output is useful to make the frame intervals as short as possible in external trigger mode. For example, when the camera is working in the edge trigger mode, the next frame can start after the previous frame exposure is done. Thus, the camera cannot accept a trigger for the next frame during the exposure period. To reduce useless time to be as short as possible, it is necessary to know the period when the camera can accept a trigger for the next frame. The trigger ready output shows the trigger ready period when the camera can accept an external trigger in the external trigger mode.

10-2-7 SLOW SCAN

It can change from the standard scan (Default) to the slow scan with a command. Read-out noise is reduced with slow scan.

Note

- Refer to Figure 10-2-2 FRAME RATE CALCULATION about the frame rate value and calculation formula with slow scan.

10-2-7-1 Timing Diagram

The read-out time is 33 ms at full frame and the frame rate is 30 frames/s at fastest. The timing diagram for slow scan is similar to the timing diagram for the standard scan except readout time. Please refer to the timing diagram of standard scan for each trigger mode.




Figure 10-13

10-2-8 GLOBAL RESET

10-2-8-1 Timing diagram

Global reset function enables to reset the electric charge of all pixels at the same time. Then all pixels can start exposure at the same time.
Global reset can work with Edge trigger mode and Level trigger mode.

[Edge trigger mode]




Figure 10-14

[Level trigger mode]




Figure 10-15

10-3 LIGHT SHEET READOUT MODE

Lightsheet Readout Mode is a unique feature of CMOS image sensor which provides improved control over the rolling shutter mechanism.

By finely synchronizing the camera readout with the illumination scan, scattered light is rejected allowing images of higher signal to noise ratios to be acquired.

The detail information of lightsheet readout mode is published on our website.

Website

<http://www.hamamatsu.com/jp/en/technology/innovation/lightsheetreadout/index.html>

10-3-1 READ-OUT DIRECTION

With normal area mode, the camera readout from the center line to the top line and to the bottom line simultaneously.

With light sheet readout mode, the camera readout from the top to the bottom line or from the bottom to the top line.

Light sheet readout mode works with only Camera Link interface operation. And the fastest frame rate is 49 frames /s.




Figure 10-16 Normal area mode




Figure 10-17 Light sheet readout mode

- (1) Top to bottom read-out: The data is read-out from the top to the bottom line.
- (2) Bottom to top read-out: The data is read-out from the bottom to the top line.




Figure 10-18 Top to bottom read-out




Figure 10-19 Bottom to top read-out

Light sheet readout mode supports three operation modes. Also light sheet readout mode supports sub-array read-out.

- (1) Free running mode
- (2) Edge trigger mode
- (3) Start trigger mode

10-3-2 FRAME RATE CALCULATION

[Calculation formula]

$$1/(\text{Exp1} + (\text{Vn} + 10) \times 1\text{H})$$

Vn = Number of vertical line (the number of lines from the starting edge of the sensor.)

Exp1 = 9.744 36 μ s to 10 s

1H = 9.744 36 μ s to 100 ms



- The Exp1 value must be input to the calculation formula below in units of seconds.

[Value of frame rate]

Horizontal width × Vertical width		Frame rate (frame/s)
2048	2048	49
	1024	99
	512	196
	256	384
	128	738
	64	1368
	8	5401
	4	6841

Size and position for sub-array read-out can be configured according to the table below.

Settings	
Size	Vertical 4 lines steps
Position	

10-3-3 READ-OUT TIME OF THE HORIZONTAL LINE

Read-out time and exposure time can be varied with light sheet readout mode.

$$\text{Readout time} = \text{Vn} \times 1\text{H}$$

Vn = Number of vertical line

1H = 9.744 36 μ s to 100 ms

The range of exposure time with light sheet readout mode is as shown below.

The maximum exposure time can be decided according to the readout time of a frame, which is longer than 10 s or not.

$$\begin{array}{ll} 1\text{H to } (\text{Vn} \times 1\text{H}) & \text{when } (\text{Vn} \times 1\text{H}) < 10 \text{ s} \\ 1\text{H to } 10 \text{ s} & \text{when } (\text{Vn} \times 1\text{H}) \geq 10 \text{ s} \end{array}$$

10-3-4 TIMING DIAGRAM

[Free running mode]




Figure 10-20 Free running mode

[Edge trigger mode]




Figure 10-21 (Ex. rising edge)

[Start trigger mode]




Figure 10-22 (Ex. rising edge)

10-3-5 PROGRAMMABLE TIMING OUTPUT

By using the programmable timing output, synchronizing with external devices is simple. A system which needs simple timing signal does not require a delay unit or pulse generator. It is possible to program and output a pulse that has an optional pulse width and an optional delay time to the end of readout timing, Vsync or Hsync. The range of delay is 0 μ s to 10 s, and the range of pulse width is 1 μ s to 10 s.

The relation between the parameter which can be set with each reference signal, and an output signal becomes as shown below.

Reference signal	Output signal
Read End	Camera outputs a pulse after certain delay from the end of sensor readout. Also the pulse width can be set.
Vsync	Camera outputs a pulse after certain delay from the beginning of readout. Also the pulse width can be set.
Hsync	Camera outputs a pulse after certain delay from the ending of readout. Also the pulse width can be set.




Figure 10-23 Programmable timing output




Figure 10-24 Programmable timing output referenced with Hsync

When you choose Hsync for the reference of programmable timing output, camera can output some pulses before start the exposure. It is called as Pre-Hsync. You can set the number of Pre-Hsync with a software command.

10-4 REAL-TIME CORRECTION FUNCTIONS

There are a few pixels in CMOS image sensor that have slightly higher readout noise performance compared to surrounding pixels. The camera has real-time variant pixel correction features to improve image quality. The correction is performed in real-time without sacrificing the readout speed at all. This function can be turned ON and OFF. (Default is ON)

11. PRECAUTION WHEN USING CMOS IMAGE SENSOR

This camera uses scientific CMOS image sensor. Careful attention must be paid to the following points when using CMOS image sensor:

(1) White spot

Subjecting CMOS image sensor to extended exposures may cause failure in part of the silicon wafer, resulting in white spots. Currently this phenomenon is not currently preventable. If CMOS image sensor is at a fixed temperature, recurrence of the white spot increases proportionally with the exposure time, so this can be rectified with dark subtraction*. Cosmic ray may generate white spot.

* After acquiring an image using a certain exposure time is loaded, CMOS image sensor is exposed to darkness for the same amount of time, and another image is obtained. After this, the difference between the images is determined, and the data for the dark portion of the original image is nullified.

(2) Folding distortion

A rough-edged flicker may be visible when imaging striped patterns, lines, and similar subject matter.

(3) Over light



- Be careful not to input too strong light such as high-energy laser into CMOS image sensor because CMOS image sensor may be damaged by over light.

12. MAINTENANCE

12-1 CARE

Perform cleaning of this equipment with the dry soft cloth.



- Do not wipe with a damp cloth or unclean cloth.

Then, the glass window on the image sensor should be cleaned according to the following.

- (1) Blow the dust from the glass window with an air duster.
- (2) Moisten a lens cleaning paper with a little ethanol, and wipe over center area of the window, gently.



- Use Lens Cleaning Paper for cleaning of glass window in front of the image sensor.



- Use a plastic tweezers and take extra care not to scratch the glass window with the tweezers. Even with plastic tweezers, there is possibility to make scratch on the glass window in case tweezers touch it.



- Avoid touching the surrounding parts of image area when wiping the glass window.

- (3) Confirm whether dust is not left.

Attach the camera to an optics, and check if there is dust or not under the uniform light condition. If there is dust on the image, please clean the glass window again.

12-2 INFORMATION ON COOLING WATER FOR THE CIRCULATING WATER COOLER



- Regarding handling cooling water and circulating water cooler, refer to instruction manual attached to the circulating water cooler.



- It is recommended to use soft water (except pure water) for cooling water.



- Do not use hard water for cooling. It cause inside of cooling water circulating path to be calcified or corroded and it result lower flow rate or water flow stop. When using hard water, conduct a process to soften water before use it.

12-2-1 WHEN USING COOLING WATER OTHER THAN RECOMMENDED

Note

Pure water

- Pure water is not appropriate for cooling water. There is possibility that pure water absorb component of cooling water path and it may cause corrosion. In addition pure water is easy to be polluted and cause impurity, sliminess or forming foreign substances. It cause lower flow rate or water flow stop.

Note

Distilled water / Deionized water

- When using the camera inside clean room, it is possible to use distilled water or deionized water by conducting periodical check. However notice it increases possibility of corrosion inside cooling water path, lowering flow rate or water flow stop.
 - Monthly check : Check water impurity, non-existence of sliminess, foreign particle is not mixed with water or not adhered inside water path and no unusual odor. If you find any of the issues, exchange cooling water and clean cooling water path.

Note

Soft water from tap

- It is possible to use soft water from tap with conducting periodical change of cooling water and checkup. However notice it increases possibility of corrosion inside cooling water path, lowering flow rate or water flow stop.
 - Monthly check : Check water impurity, non-existence of sliminess, foreign particle is not mixed with water or not adhered inside water path and no unusual odor. If you find any of the issues, exchange cooling water and clean cooling water path.
 - Exchange cooling water every 3 months.
 - Clean cooling water path every 6 months.

Note

Bottled water

- One example of soft water which is commonly available is mineral water (Hardness less than 70). Check hardness of water by referring product information of bottled water manufacturer.

13. TROUBLESHOOTING CHECKLIST

If an abnormality occurs, look up the possible causes in the following tables and, if necessary, report the details to Hamamatsu subsidiary or local distributor.

13-1 IMAGE IS NOT TRANSFERRED

Cause	Measures	Chapter
AC adapter or other cable is loose	Reconnect the cable	7
AC adapter or other cable is broken	Replace the cable	7
The correct command has not been sent to the camera	Recheck command	

13-2 ALTHOUGH IMAGES ARE TRANSFERRED

(1) Scratches or discoloration visible on the screen

Cause	Measures	Chapter
Lens is dirty	Wipe the lens	12

(2) Image is blurred

Cause	Measures	Chapter
Lens is not focused	Contact Hamamatsu subsidiary or local distributor	17
Condensation appear	Confirm the operating environmental conditions	8

(3) Only shadowed images are output

Cause	Measures	Chapter
Lens mount cap has been left on	Remove the cap	
Amount of light is too much or too low	Reduce amount of light	

(4) All screens overflow

Cause	Measures	Chapter
Too much amount of light	Reduce amount of light	
Contrast enhancement is too high	Reduce gain	

(5) Noise appears on the screen

Cause	Measures	Chapter
Exogenous noise	Find and remove cause	
Poor connection of internal connector	Contact Hamamatsu subsidiary or local distributor	17
Defective circuit system		

14. SPECIFICATIONS

14-1 CAMERA SPECIFICATIONS

(1) Electric specifications

		Camera Link	USB 3.0
Imaging device		Scientific CMOS image sensor	
Effective number of pixels		2048 (H) × 2048 (V)	
Cell size		6.5 μm (H) × 6.5 μm (V)	
Effective area		13.312 mm (H) × 13.312 mm (V)	
Full well capacity (typ.)		30 000 electrons	
Readout noise (typ.)		1.0 electrons (median) 1.6 electrons (r.m.s) at standard scan 0.8 electrons (median) 1.4 electrons (r.m.s) at slow scan	
Dynamic range *1		37 000 : 1	
Cooling method		Peltier device + Forced air-cooled, Water-cooled	
Cooling temperature	at Forced air-cooled	- 10 °C (Ambient temperature: + 20 °C)	
	at Water-cooled	- 20 °C (Water temperature: + 20 °C)	
	at Maximum cooling	- 30 °C (typ.) (at Water-cooled; Water temp. 15 °C, Ambient temp. 20 °C)	
Frame rate Top: standard scan Bottom: slow scan	at Full resolution	100 frame/s 30 frame/s	30 frame/s
	at 1024 lines at center position	200 frame/s 60 frame/s	60 frame/s
	at 8 lines at center position	25 655 frame/s 7 696 frame/s	7894 frame/s
	at Horizontal 512 pixels at 8 lines at center position	---	25 655 frame/s
Readout mode	Normal readout mode	1×1	
	Binning readout mode	2×2,4×4 (Digital binning) *2	
	Sub-array readout mode	Upper and lower symmetric area	Configurable for each vertical 8 pixels and horizontal 512 pixels.
Light sheet readout mode	Readout time	20 ms to 204.8 s	---
	Exposure time	9.7 μs to 10 s	---
	Readout mode	Full/Sub-array	---
Dark current	at - 10 °C	0.06 electron/pixel/s	
	at - 20 °C	0.02 electron/pixel/s	
	at - 30 °C	0.006 electron/pixel/s	
A/D bit		16 bit	
Linearity		Less than 3 %	
Exposure time (Standard scan)	Free running mode	1 ms to 10 s	
	Free running mode / Sub-array mode	38.96 μs to 10 s	
	External trigger mode	1 ms to 10 s	

		Camera Link	USB 3.0
Exposure time (Slow scan)	Free running mode	3 ms to 10 s	
	Free running mode / Sub-array mode	129.99 μ s to 10 s	
	External trigger mode	3 ms to 10 s	
Conversion factor		0.46 electrons/count (typ.)	
Dark offset		100 counts (at Normal readout)	
External trigger Function		Edge trigger, Level trigger, Synchronous readout trigger, Start trigger	
External signal input		External input (SMA connector) (Camera Link)	External input (SMA connector)
External trigger input level		3.3 V LVCMOS level	
External trigger output level		3.3 V LVCMOS level	
External trigger delay function		0 μ s to 10 s (10 μ s steps)	
External signal output		Global exposure timing output, Trigger ready output, Programmable timing output 1, Programmable timing output 2, Programmable timing output 3 (Continuous High or Low output)	
Image processing function		Real-time offset correction, Real-time gain correction, Real-time defect pixel correction (Default)	
Interface		Camera Link Full configuration ^{*3}	USB 3.0 Super Speed
Lens mount		C-mount (C11440-22CU) / F-mount (C11440-22CU01)	

* 1 Calculated from the ratio of the full well capacity and the readout noise.

* 2 Digital binning processing in the camera.

* 3 Original mode based on 80 bit mode

(2) Power supply specifications

Camera	Input power supply	DC 12 V
	Typical output	---
	Power consumption	55 W
AC adapter	Input power supply	AC 100 V to AC 240 V 50 Hz / 60 Hz 2.5 A
	Typical output	DC 12 V 8.34 A
	Power consumption	70 VA

Note

- Fluctuations of input power supply voltages are not to exceed $\pm 10\%$ of the nominal voltage.

(3) Operating environment

Ambient operating temperature	0 °C to + 40 °C	
Ambient storage temperature	-10 °C to + 50 °C	
Ambient operating humidity	at Forced air-cooled	Less than 70 %, no condensation
	at Water-cooled	Less than 50 % (at ambient temperature 20 °C)
	at Maximum cooling	Less than 40 % (at ambient temperature 20 °C)
Ambient storage humidity	Less than 90 %, no condensation	
Operating space	Indoor, altitude up to 2000 m	

(4) Dimensional outline and weight

Dimensional outline	85 mm (W) × 85.5 mm (H) × 125 mm (D)
Weight	Approx. 2.0 kg (Camera only)

Note

- See Chapter 15 [DIMENSIONAL OUTLINES] for detail of dimensions.



- Be careful not to drop off the camera or not drop underfoot when making it move because it is approx. 2.0 kg.

(5) Applicable standards

EMC	EN61326-1: 2013 Class A
-----	-------------------------

14-2 CONDENSATION

At the Water-cooling, if ambient temperature and ambient humidity become high, condensation will take place easily. Use the camera under the environment where condensation will not take place referring to the following graph.




Figure 14-1

14-3 SPECTRAL RESPONSE CHARACTERISTICS(TYP.)




Figure 14-2

14-4 INTERFACE SPECIFICATIONS

14-4-1 CAMERA LINK INTERFACE

The camera is based on Camera Link interface [Camera Link full configuration Deca mode], that transfers 80 bit data in parallel (8 bit x 10 port). It is an extended interface of Camera Link full configuration Standard.

(1) Pin assignments of Camera Link connector 1 (SDR-26)

Camera connector	Frame grabber connector	Channel Link signal
1	1	Inner Shield
2	25	X0-
3	24	X1-
4	23	X2-
5	22	Xclk-
6	21	X3-
7	20	SerTC+
8	19	SerTFG-
9	18	CC1-
10	17	CC2+
11	16	CC3-
12	15	CC4+
13	13	Inner Shield
14	14	Inner Shield
15	12	X0+
16	11	X1+
17	10	X2+
18	9	Xclk+
19	8	X3+
20	7	SerTC-
21	6	SerTFG+
22	5	CC1+
23	4	CC2-
24	3	CC3+
25	2	CC4-
26	26	Inner Shield

(2) Pin assignments of Camera Link connector 2 (SDR-26)

Camera connector	Frame grabber connector	Channel Link signal
1	1	Inner Shield
2	25	Y0-
3	24	Y1-
4	23	Y2-
5	22	Yclk-
6	21	Y3-
7	20	Terminated
8	19	Z0-
9	18	Z1-
10	17	Z2-
11	16	Zclk-
12	15	Z3-
13	13	Inner Shield
14	14	Inner Shield
15	12	Y0+
16	11	Y1+
17	10	Y2+
18	9	Yclk+
19	8	Y3+
20	7	Terminated
21	6	Z0+
22	5	Z1+
23	4	Z2+
24	3	Zclk+
25	2	Z3+
26	26	Inner Shield


(3) Camera Link bit assignments

28 bit solution pin name	Port	Plug No.1, Channel Link X	Port	Plug No.2, Channel Link Y	Port	Plug No.3, Channel Link Z
TxIN0	Port A0	D0_0	Port D2	D1_10	Port G5	D3_5
TxIN1	Port A1	D0_1	Port D3	D1_11	Port G6	D3_6
TxIN2	Port A2	D0_2	Port D4	D1_12	Port G7	D3_7
TxIN3	Port A3	D0_3	Port D5	D1_13	Port H0	D3_8
TxIN4	Port A4	D0_4	Port D6	D1_14	Port H1	D3_9
TxIN5	Port A5	D0_5	Port D7	D1_15 (MSB)	Port H2	D3_10
TxIN6	Port A6	D0_6	Port E0	D2_0	Port H3	D3_11
TxIN7	Port A7	D0_7	Port E1	D2_1	Port H4	D3_12
TxIN8	Port B0	D0_8	Port E2	D2_2	Port H5	D3_13
TxIN9	Port B1	D0_9	Port E3	D2_3	Port H6	D3_14
TxIN10	Port B2	D0_10	Port E4	D2_4	Port H7	D3_15 (MSB)
TxIN11	Port B3	D0_11	Port E5	D2_5	Port I0	D4_0
TxIN12	Port B4	D0_12	Port E6	D2_6	Port I1	D4_1
TxIN13	Port B5	D0_13	Port E7	D2_7	Port I2	D4_2
TxIN14	Port B6	D0_14	Port F0	D2_8	Port I3	D4_3
TxIN15	Port B7	D0_15 (MSB)	Port F1	D2_9	Port I4	D4_4
TxIN16	Port C0	D1_0	Port F2	D2_10	Port I5	D4_5
TxIN17	Port C1	D1_1	Port F3	D2_11	Port I6	D4_6
TxIN18	Port C2	D1_2	Port F4	D2_12	Port I7	D4_7
TxIN19	Port C3	D1_3	Port F5	D2_13	Port J0	D4_8
TxIN20	Port C4	D1_4	Port F6	D2_14	Port J1	D4_9
TxIN21	Port C5	D1_5	Port F7	D2_15 (MSB)	Port J2	D4_10
TxIN22	Port C6	D1_6	Port G0	D3_0	Port J3	D4_11
TxIN23	Port C7	D1_7	Port G1	D3_1	Port J4	D4_12
TxIN24	LVAL	LVAL	Port G2	D3_2	Port J5	D4_13
TxIN25	FVAL	FVAL	Port G3	D3_3	Port J6	D4_14
TxIN26	Port D0	D1_8	Port G4	D3_4	Port J7	D4_15 (MSB)
TxIN27	Port D1	D1_9	LVAL	LVAL	LVAL	LVAL
TxCLKIn	PClk	Pixel Clock A,B,C	PClk	Pixel Clock D,E,F	PClk	Pixel Clock G,H,I,J


LVAL (Line Valid signal)	This signal show the period during which the line part of the image data from CMOS image sensor is in effect. This is "ON" when during the period the line is active.
FVAL (Frame Valid signal)	This signal shows the period during which the vertical part of the image data from CMOS image sensor is in effect. This is "ON" during the period the frame is active.
D0_0 to D4_15 (Digital image data)	This is the image signal data from CMOS image sensor converted A/D. "D0 to D4" has 16 bit data in each. MSB shows the most significant bit. Please see 14-4-2 [OUTPUT TIMING SPECIFICATIONS] for details.

14-4-2 OUTPUT TIMING SPECIFICATIONS


(1) Normal readout




(2) 2x2 binning readout



(3) 4x4 binning readout




(4) Sub-array readout




15. DIMENSIONAL OUTLINES

15-1 C11440-22CU



15-2 C11440-22CU01



16. WARRANTY

Hamamatsu Photonics have fully inspected this camera and checked that its performance conforms to specifications. In the unlikely event of breakdown or other malfunction, contact Hamamatsu subsidiary or local distributor.

- (1) Unless otherwise stated by Hamamatsu subsidiary or local distributor, this system is under warranty for 24 months from the delivery date.
 - Degradation with cosmic rays, the radiation (X-rays, gamma rays, UV light, etc.) of CMOS image sensor is excepted.
- (2) The warranty only covers defects in the materials and manufacturing of the system. You may be liable for repairs during the warranty period in the event of a natural disaster or if you handle the system contrary to the instructions in this manual, use it without due caution, or try to modify it.
- (3) We will repair the system or replace it, subject to availability, free of charge within the terms of the warranty.

REPAIRS

- (1) If you notice anything wrong with the camera, confirm whether or not it is malfunctioning by referring to the troubleshooting checklist in this instruction manual. You must first clarify the symptoms in order to avoid any misunderstanding or error.
- (2) If you have any trouble or are unclear about anything, contact Hamamatsu subsidiary or local distributor giving the product name, serial number and details of the problem. If Hamamatsu Photonics consider the problem to be a malfunction, we will decide whether dispatch an engineer or have the camera returned to us for repairs.

17. CONTACT INFORMATION

Manufacturer

HAMAMATSU PHOTONICS K. K., Systems Division

812 Joko-cho, Higashi-ku, Hamamatsu City, 431-3196, Japan

Telephone (81) 53-431-0124, Fax: (81) 53-435-1574

E-mail: export@sys.hpk.co.jp

Local contact information worldwide could be found under:

www.hamamatsu.com

- *The contents of this manual are subject to change without notice.*
- *The unauthorized reproduction or distribution of parts or all of this manual is prohibited.*
- *If one of the following problems occurs, please contact Hamamatsu Photonics.
(See the CONTACT INFORMATION.) We will deal with the problem immediately.*
 - *Some contents of the manual are dubious, incorrect or missing.*
 - *Some pages of the manual are missing or in the wrong order.*
 - *The manual is missing or dirty.*